


Searching via Source Search



The screenshot shows the JAGGAER Source Search interface. The 'Source Search' tile is highlighted with a red box. A search for 'acetone' is shown with filters for 'Internal', 'External', and 'All'. The 'Labs' source channel is selected and highlighted with a red box. Below, the search results for 'Acetone-2,4-dinitrophenylhydrazone' and 'Acetone' are displayed, with the 'Labs' source channel selected for the first result.

Search Term: acetone

Filters:

- Labs
- Stockrooms
- Storerooms
- Hosted Catalogs
- eMolecules

Results:


Source Channel	Item Name	Bar Code	CAS #	MDL #	Quantity
Labs	Acetone-2,4-dinitrophenylhydrazone	19761	1567-89-1	MFCD00137172	3 G (3 G)
Labs	Acetone	346043	67-56-1		

- ▶ Select the “Source Search” tile from the home screen
- ▶ Select the “Labs” option for searching
- ▶ Enter search parameter into the “search term” text field
- ▶ Click the “Search” button or press the “enter” key
- ▶ Select the “Labs” source channel to view the results.

Note: From the search result window, user can add items from “Labs” source to a “To Do”.

Searching laboratory Inventory

Searching via Lab Inventory



JAGGAER ERM Researcher

Menu | Home Page

Source Search | Structure Search | Shopping Cart

To Do List | My Requests | Lab Inventory

Receiving | Container Operations | Material Search

Create Material | Held Carts | Open Approval Items

JAGGAER ERM Researcher

Menu | Home Page | Containers and Materials | Lab Inventory

Owner: OR Location: CHEMISTRY>F10>3>333

Search

	Material	Bar Code	CAS #	Package	Status
<input checked="" type="checkbox"/>	'L-beta-Homoprolin' hydrochloride	74897	53912-85-9	250.0000153 MG	AVAILA
<input type="checkbox"/>	(+)-2,2'-Isopropylidenebis[(4R)-4-phenyl-2-oxazolin]	151538	150529-93-4	1 g	INTRAN
<input type="checkbox"/>	(+)-8-Methoxydiisopinocampheylborane	152803	99438-28-5	5 g	AVAILA

- ▶ Select “Lab Inventory” tile.
- ▶ Select the location in the drop-down menu
- ▶ Click “Search” button, which brings up the container in the locations
- ▶ To edit or print barcode label for container, click on the “Edit Container” icon

Note: The drop-down location list is link to the user’s “Alternate Delivery Location”.